


No.1 IN THE WORLD IN TERMS OF EFFICIENCY

Vaginal Veil Collector V-Veil UP2™



Patents pending


Comparison of effectiveness between two methods




Almost twice the ability to detect high-risk oncogenic papillomaviruses (HPV-HR) by multiplex real-time PCR, using the V-Veil UP2 Vaginal Veil Collector, compared to the standard swab sampling performed by a gynaecologist or a nurse.


— Instructions for use of V-Veil UP2 —

Follow the instructions below before using this medical device::

- Wash your hands properly before using this medical device
 - Do not use if you have had sexual activities more recently than 2 days
 - Do not use if you have washed the vaginal cavity more recently than 48 hours
 - If you are pregnant or you have blood in the vaginal cavity use this medical device only at the recommendation of your doctor.
- 



No restrictions during these 2 to 30 minutes before removing.



Vaginal Veil Collector

Patent pending



V-shaped
for
collect




Advantage of flaps collector: our special flaps design collect and retain more cells from the cervix area.



STUDY

The graphic below compares the molecular detection of HPV genotypes performed from:

- Genital secretions self-collected by women themselves, using the Vaginal Veil Collector V-Veil Up2
 Genital secretions collected by a gynaecologist or a nurse, using a swab



This study was performed on 256 adult women by the research team of the Laboratory of Virology of the Hospital Georges Pompidou in France. The veil detects 12 times more oncogenic HPV-HR than the use of a standard swab performed by a gynaecologist. (Details on study)

No.1 IN THE WORLD IN TERMS OF EFFICIENCY

Vaginal Veil
Collector

حجاب مهپلی جامع

Copyright V-Veil-UP

V-Veil UP2

14 mm Veil

14 mm



Box of 1 unit

Adapted for Ø 15 mm storage container

Box of 250 units



x1

x1



+



UV-C Disinfected



x250

x1



+



UV-C Disinfected

V-Veil UP2

28 mm Veil

28 mm



Box of 1 unit

Adapted for Ø 30 mm storage container

Box of 200 units



x1


x1



+



UV-C Disinfected



x200

x1



+



UV-C Disinfected

UP2 KIT

28 mm Veil

28 mm



Box of 1 KIT UP2



x1

x1

x1



+



+



UV-C Disinfected UV-C Disinfected

Adapted for Ø 15mm storage container
thanks to the 2nd applicator
serving as a pusher

Produced in an EN ISO 13485:2016 certified facility. Registered at European database on medical devices (EUDAMED)
Authorized by National Agency for Medicines and Medical Devices of Romania